CELLBANKER® 2

Cryopreservation Medium (Serum Free)

Cat # 11891 Qty: 100ml

Expiry Date: 3 years from manufacturing date (see label)




Cell Freezing:

For optimum results, cells for cryopreservation should be in log phase of growth. Similar or standard freezing protocols may be substituted.

- 1. Examine and make sure the cell culture is free of contamination, in healthy situation and proper confluency, etc.
- 2. Perform a cell count to determine the viability of cells
- 3. Gently pellet the cells by centrifugation (3 5 minutes at 1,000~2,000rpm, 4°C). Remove the supernatant by using an aspirator.
- 4. Gently suspend the cells with CELLBANKER[®] 2 cryopreservation medium (1 ml for $5 \times 10^5 5 \times 10^6$ cells).
- 5. Dispense the cell suspension in 1ml aliquots to cryopreservation vials that have been labeled with the cell line name, cell concentration, passage date and other essential information.
- 6. Place the vials directly in a -80°C for storage. If necessary, transfer the frozen vials to a liquid nitrogen storage tank after the vials have been frozen for at least 24 hours.
- 7. Optimum protocol may change with the cell types.

IMPORTANT: Optimum protocol may change with the cell types.

Procedure for Use:



Thawing:

- 1. Remove the frozen cell from storage and quickly thaw in a 37°C shaking water bath.
- 2. Immediately dilute and gently mix each 1ml of cells with 10ml of complete cell culture medium.
- 3. Gently pellet the cells by centrifugation (3-5 minutes at 1,000 2,000rpm, 4°C). Remove the supernatant by aspirator.
- 4. Gently suspend the cells with appropriate volume of complete cell culture medium and plate in a culture flask.
- 5. Continue the further culture procedures according to standard protocols.

Guarantee of Quality:

- 1. Bacterial contamination free Product has been tested and confirmed to be free of bacteria, fungi and mycoplasma.
- 2. Chemical Analysis: pH (7.0 to 8.5 at room temperature) Endotoxin (<5 EU/mL)
- 3. Performance test Cell viability above 80% (JM404, SK-007) is guaranteed.

Storage of CELLBANKER® 2:

- 1. CELLBANKER® 2 should be stored at 4°C or below.
- 2. Keep frozen in case of long term unused period (3 months or longer)
- 3. Repeated freezing and thawing may impair the quality of the product; it is recommended that CELLBANKER® 2 is aliquoted before freezing.

Precautions:

- 1. For research use only
- 2. Not for clinical or diagnostic use.
- 3. Performance of trial tests using cells of intended use before experiments is recommended.

References:

Bagley. J. A. et al. (2017) Fused cerebral organoids model interactions between brain regions. Nature Methods, doi:10.1038/nmeth.4304

Product Range:

Description	Pack Size
CELLBANKER® 1 - Serum Containing	20 ml
CELLBANKER® 1 - Serum Containing	4 x 20 ml
CELLBANKER® 1 - Serum Containing	100 ml
CELLBANKER® 2 - Serum Free	20 ml
CELLBANKER® 2 - Serum Free	4 x 20 ml
CELLBANKER® 2 - Serum Free	100 ml
STEM-CELLBANKER® - GMP	20 ml
STEM-CELLBANKER® - GMP	4 x 20 ml
STEM-CELLBANKER® - GMP	100 ml
STEM-CELLBANKER® - GMP - DMSO Free	20 ml
STEM-CELLBANKER® - GMP - DMSO Free	4 x 20 ml
STEM-CELLBANKER® - GMP - DMSO Free	100 ml
CELLOTION cell wash solution	100 ml

Cell Types Tested:

Cell type	Description				
293	Sheared human Ad5 DNA-transformed cell line.				
293T	Human cell line expressing SV40 large T antigen				
32D	Murine myeloid cell				
3LL	Murine Lewis lung cancer cell				
A10	Thoracic aortic smooth muscle of an embryonic rat				
A2781	Glioblastoma cell				
Ac2F	Rat liver cell				
AtT20	1 ******				
Ba/F3	Murine anterior pituitary cell Murine Pro-B cell				
BHK-21	1				
C2C12	Baby Hamster renal fibroblast cell				
CZC1Z	Mouse myoblast cell				
C3H10T1/2	C3H10T1/2 Murine embryonic fibroblast cell				
Caco-2	Human colonic adenocarcinoma cell				
CHL/1U	Neonatal Chinese Hamster Lung cell				
СНО	Chinese Hamster Ovary cell				
CHO-K1	Chinese Hamster Ovary cell				
Colo203	Human colonic adenocarcinoma cell				
COS1	African green monkey kidney cell				
COS7	African green monkey kidney cell				
CTLL-2	Murine T-cell				
DLP-1	Human colonic adenocarcinoma cell				
DT40	Chicken B lymphocyte cell				
DU145	Human prostate carcinoma cell				
EBV transf	ormed B cell				
EJ-1	Human bladder carcinoma				
ELM-D	Murine erythroblastic leukemia cell				
HeLa cell	Human Uterine Cervical Carcinoma Cell				
Нер3В	Human hepatoma cell				
Hepal-6					
Hepatocyt	·				
HepG2	Hepatocellular carcinoma cell				
HL-60	Human promyelocytic leukemia cell				
HT-2	Murine cell				
Huh-7					
Jurkat	Human hepatocarcinoma cell				
K-562	Human T-cell line				
Kato-III	Human Caucasian Chronic myelogenous leukaemia cell Human stomach carcinoma				
Keratinocy					
KM12-LX	Human Caucasian Chronic myelogenous leukaemia cell				
L5178Y	Murine lymphoma cell				
L929	Murine fibroblast cell				
LM	Murine fibroblast cell				
LNCap	Human prostate adenocarcinoma cell				
Lymphocy					
MCF-7	Human metastatic mammary carcinoma cell				
MDCK					
Microvascular Endothelial Cell					
Molt-4	Human T-cell				
NB- 1/GOTO	Human neuroblastoma				
NCI-H441	Human lung adenocarcinoma epithelial cell				
NIH3T3	Mouse embryonic fibroblast cell				

P388	Murine leukimia cell			
P3/x63- Ag8.U1	Murine myeloma cell			
PC12	Rat pheochromocytoma cell			
Periodontal ligament fibroblast cell				
Periodontal ligament membrane fibroblast cell				
Raji	Human B-cell line			
RAW264.7	Murine Macrophage cell			
sf9	Insect cell			
SK-N-MC	Neuroepithelioma cell			
SN12C	Renal carcinoma cell. Subclone: PM6,Clone8, MM3			
Stomach carcinoma cell				
THle3	Human liver cell			
T24	Human bladder carcinoma			
Vero	African green monkey kidney epithelial cell			
WEHI3B	Murine myeloid leukemia			
WiDr	Human colon adenocarcinoma cell			

Typical Experimental Results:

CUT	5	Viability of cells (%)		
Cell Type	Preservation period (year)	-80°C	-196°C	
Mouse	·			
Hybridoma	10	95	95	
Myeloma	10	90	90	
L929	10	90	90	
FM3A	5	90	90	
BALB/3T3	5	90	90	
M1	5	90	90	
YAC-1	5	90	-	
Rat				
RLC-16	5	90	90	
NRK	5	90	90	
PC-12	5	90	-	
Hamster				
CHO	5	90	90	
V79	5	90	90	
Monkey				
COS-1	5	90	90	
Vero	5	90	90	
Human				
Kidney-derived tumor cell	5	90	90	
EBV transformed cell	5	90	90	
HEL-derived fibroblast	5	90	90	
Melanoma	5	90	90	
Caco-2	3	90	-	
C-5	5	90	90	
CEM	5	90	90	
K562	10	90	90	
Jurkat	10	90	90	
BALL-1	5	90	90	
HUC-Fm	5	80	80	